

Identification of free amino acids in several crude extracts of two legumes using Thin Layer Chromatography

Authors

Taghread Hudaib

Sarah Brown

Daniel Wilson

Paul E. Eady

Key words

Amino acids; Aqueous-ethanolic extract; Thin layer chromatography; Organic solvents; *Vigna unguiculata*; *Phaseolus vulgaris*.

1. Introduction

Amino acids are primary metabolites and nutritional organic compounds that are the building blocks of proteins [1]. The standard (proteinogenic) amino acids consist of an amine group (NH₂), a carboxylic acid (COOH), alpha hydrogen and a side chain [2]. Thin layer chromatography (TLC) is a simple, qualitative, sensitive and widely used method for the separation and identification of amino acids in plant extracts. Amino acids vary in their solubility in water and organic solvents depending on the nature of their side chains; water soluble amino acids have polar side chains [3].

Taghread Hudaib *
University of Lincoln, School of Life Sciences, Brayford Pool, Lincoln, LN6 7TS, UK.
Telephone number +44(0)1522 83 5326
Email: thudaib@lincoln.ac.uk
Fax: +44 1522 886974

Sarah Brown
Applied Materials Technology Ltd, Lyndon Business Park, Farrier Road, Lincoln, LN6 3RU, UK.
Email: sarah.brown@appmat.co.uk

Daniel Wilson
University of Lincoln. School of pharmacy, Brayford Pool, Lincoln, LN6 7TS
Telephone number +44 (0)1522 886825
Email: dwilson@lincoln.ac.uk

Paul E. Eady
University of Lincoln, School of Life Sciences, Brayford Pool, Lincoln, LN6 7TS, UK.
Telephone number +44(0)1522 83 5461
Email: peady@lincoln.ac.uk

In general, the 70% aqueous ethanol is the most preferable solvent to prepare crude extracts prior to amino acid TLC analysis [4].

In this study, the amino acid profiles of seed cotyledons and seed coats from two types of legumes; black-eyed beans (*Vigna unguiculata*) and red kidney beans (*Phaseolus vulgaris*) were determined and compared in order to evaluate their nutritional value as major source of plant amino acids. To improve the separation and detection of free amino acids, the TLC technique was applied on several crude extracts obtained at the primary extraction stage using a range of polar and nonpolar solvents. The obtained crude extracts are:

1. Kidney bean seeds 70% Ethanol (KBS 70% Ethanol)
2. Kidney bean seeds 100% Methanol (KBS100% Methanol)
3. Kidney bean seeds 100% Acetone (KBS 100% Acetone)
4. Kidney bean seed coat 70% Ethanol (KBSC 70% Ethanol)
5. Kidney bean seed coat 100% Methanol (KBSC 100% Methanol)
6. Black-eyed bean seeds 70% Ethanol (BEBS 70% Ethanol)
7. Black-eyed bean seeds 100% Methanol (BEBS100% Methanol)
8. Black-eyed bean seeds 100% Acetone (BEBS 100% Acetone)
9. Black-eyed bean seed coat 70% Ethanol (BEBSC 70% Ethanol)
10. Black-eyed bean seed coat 100% Methanol (BEBSC100% Methanol)
11. Kidney bean seeds 100% Chloroform (KBS 100% Chloroform)

Taghreed Hudaib *
University of Lincoln, School of Life Sciences, Brayford Pool, Lincoln, LN6 7TS, UK.
Telephone number +44(0)1522 83 5326
Email: thudaib@lincoln.ac.uk
Fax: +44 1522 886974

Sarah Brown
Applied Materials Technology Ltd, Lyndon Business Park, Farrier Road, Lincoln, LN6 3RU, UK.
Email: sarah.brown@appmat.co.uk

Daniel Wilson
University of Lincoln. School of pharmacy, Brayford Pool, Lincoln, LN6 7TS
Telephone number +44 (0)1522 886825
Email: dwilson@lincoln.ac.uk

Paul E. Eady
University of Lincoln, School of Life Sciences, Brayford Pool, Lincoln, LN6 7TS, UK.
Telephone number +44(0)1522 83 5461
Email: peady@lincoln.ac.uk

12. Black-eyed bean seeds 100% Chloroform (BEBS 100% Chloroform)

2. Experimental

2.1. Crude extracts preparation

Two kilograms of black-eyed beans and two kilograms of red kidney-beans were soaked in distilled water for 30 and 45 min respectively to remove the seed coat. Then, seeds and seed coats were left to dry at room temperature. After that, the seeds and seed coats were ground to flour using a coffee grinder and de-fatted using a nonpolar solvent (hexane). The ground seeds and seed coats were extracted using solvent extraction methods [5] in which solvents with different polarity were used in Soxhlet extraction [6]. For both the decorticated seeds and seed coats, 500 ml of solvent were needed for extraction of 100g of flour along the extraction process. The process continued until the resulting filtrate is colourless, which indicates that most of the soluble constituents have been extracted [7]. The following selective aqueous and organic solvents were used to prepare the crude extracts:

- i. Aqueous-ethanol (70% ethanol), Polarity = 1.0 for water and 0.65 for ethanol.
- ii. Methanol (100%), Polarity = 0.76
- iii. Acetone (100%), Polarity = 0.36
- iv. Chloroform (100%), Polarity = 0.26

Taghreed Hudaib *
University of Lincoln, School of Life Sciences, Brayford Pool, Lincoln, LN6 7TS, UK.
Telephone number +44(0)1522 83 5326
Email: thudaib@lincoln.ac.uk
Fax: +44 1522 886974

Sarah Brown
Applied Materials Technology Ltd, Lyndon Business Park, Farrier Road, Lincoln, LN6 3RU, UK.
Email: sarah.brown@appmat.co.uk

Daniel Wilson
University of Lincoln. School of pharmacy, Brayford Pool, Lincoln, LN6 7TS
Telephone number +44 (0)1522 886825
Email: dwilson@lincoln.ac.uk

Paul E. Eady
University of Lincoln, School of Life Sciences, Brayford Pool, Lincoln, LN6 7TS, UK.
Telephone number +44(0)1522 83 5461
Email: peady@lincoln.ac.uk

Extracts were then suction-filtered and evaporated to dryness under vacuum in a rotary evaporator (Rotavapor® R-210/R-215) at the boiling point of the solvent. The dry extracts were preserved at 4°C prior to analysis for amino acid [8].

2.2. Thin Layer Chromatography

The amino acid profiles of all the extracts were determined using the Thin Layer Chromatography technique (TLC). The following fresh solutions were prepared to run the TLC:

- i. Mobile phase: Butanol: Acetic acid glacial: Water, mixed in a ratio of 12: 3: 5 respectively [9]
- ii. Visualization reagent: Ninhydrin reagent: 0.2g/100ml ethanol.
- iii. Amino acid standards: The following Amino acid standards were dissolved in water to a concentration of 2% (w/v) and used for TLC spotting: Glycine, Serine, Leucine, Cysteine, Valine, Aspartic Acid, Tryptophan, Tyrosine, Threonine, Histidine, Proline, Glutamic Acid, Cystine, Arginine, Alanine, Glutamine, Isoleucine, Asparagine, Methionine, Phenylalanine, Hydroxyproline, and Lysine.

Using capillary tubes, small spots of amino acid standards and bean extract solutions were applied 2cm above one edge of silica coated TLC plates (Whatman silica gel 60 A, 20 × 20 cm, 250 µm thickness with fluorescent indicator) with a 1cm space between the spots. The spots were then labeled with a pencil and left to dry at room

Taghreed Hudaib *
University of Lincoln, School of Life Sciences, Brayford Pool, Lincoln, LN6 7TS, UK.
Telephone number +44(0)1522 83 5326
Email: thudaib@lincoln.ac.uk
Fax: +44 1522 886974

Sarah Brown
Applied Materials Technology Ltd, Lyndon Business Park, Farrier Road, Lincoln, LN6 3RU, UK.
Email: sarah.brown@appmat.co.uk

Daniel Wilson
University of Lincoln. School of pharmacy, Brayford Pool, Lincoln, LN6 7TS
Telephone number +44 (0)1522 886825
Email: dwilson@lincoln.ac.uk

Paul E. Eady
University of Lincoln, School of Life Sciences, Brayford Pool, Lincoln, LN6 7TS, UK.
Telephone number +44(0)1522 83 5461
Email: peady@lincoln.ac.uk

temperature before being located vertically in a TLC chamber containing the developing solution (the mobile phase). The plate edge beneath the spots was immersed in the developing solution and the TLC chamber sealed. Then, the plates were left for approximately three hours. The highest point reached by the mobile phase was marked. When the plates were completely dry they were sprayed with the Ninhydrin reagent and left at 80-100°C for 5min to allow the visualizing reagent to react and stain the spots. The distances from the starting point to the spot centers were measured and the Retardation Factor (R_f) value for each standard and unknown spot was calculated (R_f = the distance traveled by the spot / the distance traveled by the mobile phase front). The following symbols were used to describe the results of the detected amino acids: (+) observed; (++) bright; (\pm) faint [10].

3. Results and Discussion

The organic solvents used to produce the crude extracts varied in their ability to extract amino acids; no amino acids were detected in Chloroform fractions. Sixteen amino acids were extracted by different polar solvents. The calculated R_f value for each standard is shown in Table 1. These R_f values were used as markers to identify the amino acids in each extract (Table 2).

Taghreed Hudaib *
University of Lincoln, School of Life Sciences, Brayford Pool, Lincoln, LN6 7TS, UK.
Telephone number +44(0)1522 83 5326
Email: thudaib@lincoln.ac.uk
Fax: +44 1522 886974

Sarah Brown
Applied Materials Technology Ltd, Lyndon Business Park, Farrier Road, Lincoln, LN6 3RU, UK.
Email: sarah.brown@appmat.co.uk

Daniel Wilson
University of Lincoln. School of pharmacy, Brayford Pool, Lincoln, LN6 7TS
Telephone number +44 (0)1522 886825
Email: dwilson@lincoln.ac.uk

Paul E. Eady
University of Lincoln, School of Life Sciences, Brayford Pool, Lincoln, LN6 7TS, UK.
Telephone number +44(0)1522 83 5461
Email: peady@lincoln.ac.uk

100

101

104

Taghreed Hudaib *
University of Lincoln, School of Life Sciences, Brayford Pool, Lincoln, LN6 7TS, UK.
Telephone number +44(0)1522 83 5326
Email: thudaib@lincoln.ac.uk
Fax: +44 1522 886974

Sarah Brown
Applied Materials Technology Ltd, Lyndon Business Park, Farrier Road, Lincoln, LN6 3RU, UK.
Email: sarah.brown@appmat.co.uk

Daniel Wilson
University of Lincoln. School of pharmacy, Brayford Pool, Lincoln, LN6 7TS
Telephone number +44 (0)1522 886825
Email: dwilson@lincoln.ac.uk

Paul E. Eady
University of Lincoln, School of Life Sciences, Brayford Pool, Lincoln, LN6 7TS, UK.
Telephone number +44(0)1522 83 5461
Email: peady@lincoln.ac.uk

Extracts▶ Amino acids▼	KBS 70% Ethanol	KBS 100% Methanol	KBS 100% Acetone	KBSC 70% Ethanol	KBSC 100% Methanol	BEBS 70% Ethanol	BEBS 100% Methanol	BEBS 100% Acetone	BEBS 70% Ethanol	BEBS 100% Methanol
Ala		+++	+				+	±	±	+
Arg	+		+							
Asn	±									
Asp			±			+				
Cys				+					+	
Glu			+					±	±	+
Gly	+	+								
Pro			+				+			
Tyr	±	+	+				++	++	++	+
His	±	+	±							
Ile	+	+	+			+				
Leu		+++				+	±			
Lys	+	+	±			+	±			
Phe		+		+				±		
Thr	+	±		+			+			
Val	+	+++	+			+	+	±	+	+

105

106

107 The results obtained from the TLC (Table 1 and Table 2) indicate that the amino acid
108 profile of the kidney-beans and the black-eyed beans varied with different amino acids
109 detected in different crude extracts. Some amino acids such as Ala, Asp, Glu, Pro, Leu
110 and Val were poorly or not detected in the crude extracts prepared with the 70%
111 aqueous ethanol. All of the essential amino acids with the exception of Met were found
112 in the kidney-bean extracts whilst Arg, Asn, Met, Gly and His were absent from the
113 black-eyed bean extracts. Most of the amino acids were detected in the seed

Taghreed Hudaib *
University of Lincoln, School of Life Sciences, Brayford Pool, Lincoln, LN6 7TS, UK.
Telephone number +44(0)1522 83 5326
Email: thudaib@lincoln.ac.uk
Fax: +44 1522 886974

Sarah Brown
Applied Materials Technology Ltd, Lyndon Business Park, Farrier Road, Lincoln, LN6 3RU, UK.
Email: sarah.brown@appmat.co.uk

Daniel Wilson
University of Lincoln, School of pharmacy, Brayford Pool, Lincoln, LN6 7TS
Telephone number +44 (0)1522 886825
Email: dwilson@lincoln.ac.uk

Paul E. Eady
University of Lincoln, School of Life Sciences, Brayford Pool, Lincoln, LN6 7TS, UK.
Telephone number +44(0)1522 83 5461
Email: peady@lincoln.ac.uk

cotyledon fractions, with only a few detected in the seed coat extracts. More amino acids were detected in the black-eyed bean seed coat than the kidney bean seed coat; the 70% Ethanol fraction of black-eyed bean seed was found to contain Ala, Cys, Glu, Tyr, Val and the 100% Methanol fraction contains Ala, Glu, Tyr and Val. Only one fraction of the kidney bean seed coat was found to have amino acids; the 70% Ethanol that found to contain Cys, Phe, Thr. Amino acids have different functions in different legume seeds; Cys is a sulfur containing amino acid that plays an important role in plant defense mechanisms against predators and is found to be toxic to cowpea beetles *Callosobruchus maculatus* [11]. This may explain why seed coats help protect legumes against stored seed pests and minimize losses of stored products. Moreover, both types of beans seem to lack the following amino acids: Hyp, Ser, Gln, L-Cys, Try, and Met which reduces the nutritional value of these beans. Combining beans with other legumes or grains is recommended to provide a meal that is rich in all essential amino acids.

4. Conclusion

TLC is an efficient method for the analysis of amino acids from plant materials. Amino acids vary in their polarity according to their side chains. Producing several crude extracts is a vital first step for the further separation steps. Using several solvents, with different polarity in crude extracts preparation improves the efficacy of

Taghreed Hudaib *
University of Lincoln, School of Life Sciences, Brayford Pool, Lincoln, LN6 7TS, UK.
Telephone number +44(0)1522 83 5326
Email: thudaib@lincoln.ac.uk
Fax: +44 1522 886974

Sarah Brown
Applied Materials Technology Ltd, Lyndon Business Park, Farrier Road, Lincoln, LN6 3RU, UK.
Email: sarah.brown@appmat.co.uk

Daniel Wilson
University of Lincoln. School of pharmacy, Brayford Pool, Lincoln, LN6 7TS
Telephone number +44 (0)1522 886825
Email: dwilson@lincoln.ac.uk

Paul E. Eady
University of Lincoln, School of Life Sciences, Brayford Pool, Lincoln, LN6 7TS, UK.
Telephone number +44(0)1522 83 5461
Email: peady@lincoln.ac.uk

133 the TLC technique, enhances the separation and simplifies the identification of polar
134 and non-polar amino acids.

135 **5. References:**

136 [1] *K. J. Jalkanen, M. Elstner, and S. Suhai*, J. Molecular Structure: THEOCHEM,
137 **675** (2004) 61-77.

138 [2] *A. Ambrogelly, S. Palioura, D. Söll*, Nature Chemical Biology **3** (2007) 29-35.

139 [3] *R. Wolfenden, L. Andersson, P. M. Cullis, and C. C. B. Southgate*, Biochemistry,
140 **20** (1981) 849-855.

141 [4] *B. Fried, and J. Sherma*, CRC Press, (1996).

142 [5] *J.H. Doughari*, Agricultural and Biological Sciences edition, In Tech, (2012)1-32.

143 [6] *G. Jayaprakasha, R. Singh, and K. Sakariah*, Food Chemistry, **73** (2001) 285-
144 290.

145 [7] *S.A. Audu, I. Mohammed, and H.A. Kaita*, Life Science Journal, **4** (2007) 75-79.

146 [8] *V. A. Maikai, B. V. Maikai, and P. I. Kobo*, Journal of parasitology research,
147 (2014) 904318-904318.

148 [9] *G. S. Chakraborty*, IJRAP, **1** (2010) 131-134.

149 [10] *M. Waksmundzka-Hajnos, J. Sherma, and T. Kowalska*, CRC Press, (2008).

Taghread Hudaib *
University of Lincoln, School of Life Sciences, Brayford Pool, Lincoln, LN6 7TS, UK.
Telephone number +44(0)1522 83 5326
Email: thudaib@lincoln.ac.uk
Fax: +44 1522 886974

Sarah Brown
Applied Materials Technology Ltd, Lyndon Business Park, Farrier Road, Lincoln, LN6 3RU, UK.
Email: sarah.brown@apmat.co.uk

Daniel Wilson
University of Lincoln. School of pharmacy, Brayford Pool, Lincoln, LN6 7TS
Telephone number +44 (0)1522 886825
Email: dwilson@lincoln.ac.uk

Paul E. Eady
University of Lincoln, School of Life Sciences, Brayford Pool, Lincoln, LN6 7TS, UK.
Telephone number +44(0)1522 83 5461
Email: peady@lincoln.ac.uk

150 [11] *D.H. Janzen, H.B. Juster, and E. Arthur Bell*, *Phytochemistry*, **16** (1977) 223-
151 227.
152
153
154
155
156

Taghreed Hudaib *
University of Lincoln, School of Life Sciences, Brayford Pool, Lincoln, LN6 7TS, UK.
Telephone number +44(0)1522 83 5326
Email: thudaib@lincoln.ac.uk
Fax: +44 1522 886974

Sarah Brown
Applied Materials Technology Ltd, Lyndon Business Park, Farrier Road, Lincoln, LN6 3RU, UK.
Email: sarah.brown@appmat.co.uk

Daniel Wilson
University of Lincoln. School of pharmacy, Brayford Pool, Lincoln, LN6 7TS
Telephone number +44 (0)1522 886825
Email: dwilson@lincoln.ac.uk

Paul E. Eady
University of Lincoln, School of Life Sciences, Brayford Pool, Lincoln, LN6 7TS, UK.
Telephone number +44(0)1522 83 5461
Email: peady@lincoln.ac.uk